



Dose Verification System

Patient Information Card

You have been implanted with a radiation dosimeter. This device measured the amount of radiation that you received each day during your radiation therapy. The DVS radiation dosimeter contains a ferrite core which is like a small magnet. The dosimeter is encapsulated in a hermetically sealed glass capsule.

It is important for other healthcare providers to know that you have had a dosimeter implanted. Only your healthcare provider can determine if the DVS radiation dosimeter is compatible with other therapies. If you have medical procedures performed such as an MRI (Magnetic Resonance Image), a mammogram, or other procedures, you should let your healthcare provider know about the dosimeter. Provide this card to your health care provider in advance of these medical procedures.



Dose Verification System

Patient Information Card

Direction to healthcare providers: Refer to prescribing information or contact Sichel Technologies®, Inc. for important information regarding the use of specific therapies with patients implanted with DVS® dosimeters.

Sichel Technologies®, Inc.

Contact Information:

www.dvssmartmarker.com

toll free: 1-888-387-6697/direct: 919-465-2236

email: helpdesk@sicheltech.com